

ExCell Bio

OptiVibro[®] NK Cell Serum-free Medium P01

For Research and Manufacturing Use

Not Intended for Diagnostic and Therapeutic Use

User Manual

Catalog Number NE000-N012

Catalog Number NE000-N011

Catalog Number NE000-N011S



| PRODUCT DESCRIPTION

OptiVibro® NK Cell Serum-Free Medium P01 is a serum-free, xeno -free, and cytokine-free medium specifically designed for the maintenance and expansion of natural killer (NK) cells. This medium includes OptiVibro® NK Cell Serum-free Basal Medium P01 and OptiVibro® Immune Cell Serum-free Medium Supplement UE01. Compared to traditional serum-containing media, OptiVibro® NK Cell Serum-Free Medium P01 reduces the risk of xeno-infectious agents, ensures batch-to-batch consistency, and prevents over-activation and exhaustion of NK cells due to undefined serum components, making it suitable for clinical and large-scale applications.

| SPECIFICATION, STORAGE AND TRANSPORTATION REQUIREMENT

Product Name	Cat.#	Specification	Storage	Transportation	Shelf Life
OptiVibro® NK Cell Serum-free Medium P01	NE000-N012	1000 mL kit	2-8°C Protect From Light	-	-
OptiVibro® NK Cell Serum-free Basal Medium P01	BA0092	1000 mL	2-8°C Protect From Light	< 25°C Protect From Light	12 months
OptiVibro® Immune Cell Serum-free Medium Supplement UE01	BA0332	8 mL	2-8°C Protect From Light	< 25°C Protect From Light	18 months
OptiVibro® NK Cell Serum-free Medium P01	NE000-N011	500 mL kit	2-8°C Protect From Light	-	-
OptiVibro® NK Cell Serum-free Basal Medium P01	BA0091	500 mL	2-8°C Protect From Light	< 25°C Protect From Light	12 months
OptiVibro® Immune Cell Serum-free Medium Supplement UE01	BA0331	4 mL	2-8°C Protect From Light	< 25°C Protect From Light	18 months

OptiVibro® NK Cell Serum-free Medium P01	NE000-N011S	100 mL kit	2-8°C Protect From Light	-	-
OptiVibro® NK Cell Serum-free Basal Medium P01	BA0091S	100 mL	2-8°C Protect From Light	< 25°C Protect From Light	12 months
OptiVibro® Immune Cell Serum-free Medium Supplement UE01	BA0331S	0.8 mL	2-8°C Protect From Light	< 25°C Protect From Light	18 months

| PERFORMANCE, APPLICATION AND HANDLING

RECOMMENDATIONS

1. Store cell culture medium in a dark environment, preferably in colored packaging, to protect from light exposure.
2. Avoid prolonged exposure to fluorescent or other types of lighting during transport to prevent discoloration.
3. Implement thorough cleaning and sterilization methods for transport to sterile areas; avoid UV sterilization.
4. Switch off UV lamps when transferring through UV-sterilized windows.

| INSTRUCTION FOR USE

Prepare media

1. Equilibrate OptiVibro® NK Cell Serum-free Basal Medium P01 and OptiVibro® Immune Cell Serum-Free Medium Supplement UE01 to room temperature for 1-4 hours.
2. In a biosafety cabinet, add 8 mL/4 mL of supplement to every 1 L/500 mL of basal medium. Mix by inverting 3-5 times to obtain the complete medium.

【Note】

This OptiVibro® NK Cell Serum-free Medium P01 is designed to be used in conjunction with essential cytokines to facilitate the activation and expansion of NK cells. It is crucial to supplement this medium with the necessary cytokines according to your specific cell activation protocol to achieve the desired results for NK cell

culture. The addition of cytokines provides the required signals for NK cell proliferation and function, and their absence may result in suboptimal cell growth or lack of cell activation. Always refer to your cytokine supplier's instructions for the appropriate concentrations and incubation conditions to ensure the successful expansion of NK cells.

| DISCLAIMER

1. Use the product according to the manual instructions. Deviations from these instructions are at the user's risk, and our company will not be responsible for any resulting product performance deviations.
2. This product is for scientific research and commercial production only and is not intended for clinical diagnosis or treatment. Users assume all risks for unauthorized use, and our company shall not be responsible for any consequences.